---

# **mega *LIGHTBASE***

---

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**For further Informations contact us !**

---

# mega *LIGHTBASE*

Light curing resin for fabrication of base plates  
Pink 1,4 mm

**How to use:**

The edge course for the individual base plates is marked on the plaster model. Attach space-holders if required. If no place-holders are needed, insulated the plaster model with **mega LIGHTBASE**.

Remove protective film and place the **mega LIGHTBASE** on the model, shape it and cut at the edges as marked. The transparency of the **mega LIGHTBASE** allows very accurate shaping.

Is polymerised in a normal light device (blue or UVA light). Illuminate the shaped spoon for about 5 minutes from behind.

The base plate is ready for use and further processing once it has been polymerised. Any slight inhibition layer produced by oxygen in the air can be removed by wiping the base plates with a cloth dipped in **mega Multiclean** or an alcohol-based detergent.

**Storage:**

Store in a dark place and not above 25 °C / 68°F!

**Precautions:**

In its unhardened state contains methacrylic acid ester, which can cause skin allergies. Operators who are allergic or hypersensitive should use appropriate protection. It is essential that polymerisation is sufficient and complete, as illumination for an insufficient length of time or at insufficient intensity leads to a loss of strength.

**durability:**

24 month

---

**Storage and safety measures:**

Material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapors. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

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