ERNST HINRICHSG Dental GmbH

Borsigstraße 1 · DE-38644 Goslar/Germany
Telefon +49 (0) 53 21 / 5 06 24-25 · Telefax +49 (0) 53 21 / 5 08 81
info@hinrichs-dental.de · www.hinrichs-dental.de

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DIN EN ISO 15912, Typ 2, Klasse 2

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DIN EN ISO 15912, type 2, class 2

Description:

Hinrivest® SG is a phosphate-bound, graphite-free precision investment for CrCo work. It was specially developed for the speed cast technique and is suitable for casting all types of dental alloys.

Technical data:

| | |
|--|---|
| Mixing fluid: | Hinrivest® Liquid |
| Mixing ratio, Powder : fluid: | 100 g : 18 ml |
| Mixing under vacuum: | 30 - 40 sec. |
| Working time: | 4 - 5 min |
| Storage & working temperature: | approx. 22°C (Room temperature) |
| Initial setting time: | 7 - 10 minutes |
| Setting expansion: | approx. 1.00 % (75 % Hinrivest® Liquid) |
| Preheating temperature range for placing in the furnace: | 850 - 1.000 °C (if applicable) |
| Place in the hot furnace: | 15 min after investing |

Casting models:

Weigh out the quantity required for the model and measure the Hinrivest® Liquid in compliance with the mixing ratio and required concentration. First stir by hand with a spatula and vacuum-stir for 30 - 40 seconds. Fill the surrounding parts using the stirrer at top speed. It is best to use an instrument for this process. Cast the model after mixing. Smear Vaseline over the overlapping point of the duplication band towards the side of the model to prevent cracks from forming. When duplicating models with silicone, it is unnecessary to harden and soak the models.

When duplicating models with gel, reduce the powder : liquid ratio to 100 g : 14 – 16 ml. The liquid concentration for the model and muffle should be reduced to 50 %.

Use of pressure:

Investing under pressure is **not necessary** due to Hinrivest® SG extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 15 minutes so that it does not impede the setting expansion.

Concentration of the expansion liquid: 18 ml Liquid / 100 g powder

ERNST HINRICHSG Dental GmbH

Borsigstraße 1 · DE-38644 Goslar/Germany
Tel. +49 (0) 53 21/5 06 24-25 · Fax +49 (0) 53 21/5 08 81
info@hinrichs-dental.de · www.hinrichs-dental.de

Framework with clasps: Liquid ratio Hinriwest® Liquid / distilled water: 75 % / 25 %

| 100 g | 180 g | 400 g |
|--------------|--------------|---------------|
| 14 ml / 4 ml | 24 ml / 8 ml | 54 ml / 18 ml |

Fixed/removable work: Liquid ratio Hinriwest® Liquid / distilled water: 80 % / 20 %

| 100 g | 180 g | 400 g |
|--------------|--------------|---------------|
| 14 ml / 4 ml | 26 ml / 6 ml | 58 ml / 14 ml |

Removing the duplicate model:

Carefully remove the investment model after approx. 20 min. Start waxing up once the model is cool enough to handle and the wax is no longer liable to melt.

Investing:

Measure and mix Hinriwest® SG as described above. Using the stirrer at minimum speed, pour the product until it spills from the model. Stop the stirrer, fill the muffle and stop mixing. Casting can be carried out under pressure. Mold muffles for the quick-casting procedure should be placed on a well moistened cloth when the casting is prepared under pressure. This will add a certain humidity to the already dried model.

Preheating:

Exactly 15 minutes after the casting has been poured, place the muffle with its funnel pointing towards the bottom of the furnace brought to the required pre-heating temperature (850 - 1000 °C). A maximum preheating temperature of 900 °C is recommended. This temperature is perfectly adequate for both casting the mould and final thermal expansion. The muffles will be ready for casting after preheating for 45 - 60 min, depending on their size. Keep the air circulation and exhaust systems off during the first 15 minutes and do not open the furnace (danger of bursting). Hinriwest® SG can also be preheated in a linear way at a rate of 5 - 7 °C per minute. Holding phases are no longer necessary.

Recommendations:

These technical hints are based on practical experience and experiments conducted in our R&D laboratory. The values given are purely indicative. All ERNST HINRICHSG products are subjected to strict quality controls. We reserve the right to make any technical modifications considered necessary.

Warning: Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

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